




Mr Muscle® Drain Gel

SAFE FOR USE ON ALL PIPES

Gel that clears blocked sinks, plugholes and pipes



Product Description

Sinks through standing water and blockages to clear blocked sinks, plugholes and pipes.

- 5 minute action time*
- Cuts through standing water
- Cleans pipes and controls odours
- Suitable for use in hard and / or soft water
- Safe for use on all pipes


Where to use:



Sink
(washroom)



Sink
(kitchen)



Bath

How to use:



1. Open carefully.



2. Pour the product (500ml).
Wait for 5 minutes*



3. Run drain with hot water.

DO NOT:

- Squeeze or re-use bottle
- Splash
- Plunge
- Use in toilets or sinks with waste disposal units

*Laboratory tested on bathroom clogs (hair and soap). For tougher clogs, leave longer.

Mr Muscle® Drain Gel

SAFE FOR USE ON ALL PIPES

Safety Information:

May be corrosive to metals. Causes severe skin burns and eye damage. Very toxic to aquatic life with long lasting effects. If medical advice is needed, have product container or label at hand. Keep out of reach of children. Immediately call a POISON CENTRE / doctor.

IF ON SKIN (OR HAIR): Take off immediately all contaminated clothing. Rinse skin with water / shower.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

IF INHALED: Remove person to fresh air and keep comfortable for breathing.

Store locked up. Wear protective gloves. Do not breathe fume. Dispose of contents / container in accordance with local regulations.

WARNING! Do not use together with other products. May release dangerous gases (chlorine). Keep container tightly closed in a cool place and store bottle upright.

Shelf Life

24 months from date of manufacture.

Contains:

Sodium hypochlorite, Caustic soda, Amines, C12-18-alkyldimethyl, N-oxides.

Technical Data

Appearance	Liquid
Colour	Natural
Odour	Functional
pH value	13.4 (undiluted)



Stock Code: **308557** | Case Qty: **6** | Size: **1L**